



# Test Catalog

Diagnostic. Prognostic. Predictive. Predisposition.





## NeoTYPE® Other Solid Tumor Profile

### Alternative Name

Other Solid Tumor Profile

### Methodology

Molecular

### Test Description

The NeoTYPE Other Solid Tumor Profile analyzes 30 biomarkers through a combination of next-generation sequencing (NGS), FISH, and IHC as listed below. Test orders include summary interpretation of all results to help guide treatment decisions. If Pan-TRK IHC is expressed or equivocal, reflex to either NTRK NGS Fusion Panel (Default) or NTRK 1-3 FISH Panel will be added. A microsatellite instability (MSI) NGS result of “indeterminate” will create a reflex to MSI by PCR as long as the tumor percentage is ≥20% for colorectal specimens or ≥40% with paired normal tissue available for non-colorectal specimens.

- NGS (24 genes + 2 biomarkers): AKT1, BRAF, EGFR, FGFR1, FGFR2, FGFR3, GNAS, HRAS, IDH1, IDH2, KIT, KRAS, MET, Microsatellite Instability (MSI), NOTCH1, NRAS, PDGFRA, PIK3CA, PTEN, PTPN11, SMAD4, SMO, SRC, TERT Promoter, TP53, Tumor Mutation Burden (TMB)
- FISH (2 FISH): MET, PTEN (tech-only available)
- IHC (2 biomarkers): PD-L1 LDT, Pan-TRK (tech-only available for PD-L1)

### Clinical Significance

The NeoTYPE Other Solid Tumor Profile characterizes primary or metastatic tumors of any histological subtype for the most significant genetic changes relevant to therapy decisions, prognosis, and clinical research. This test can be used for tumors arising from various tissues including liver, pancreas, prostate, kidney, head and neck, or other tumors that do not match any subtypes of the 13 other tumor-specific NeoTYPE Profiles. These other profiles are available for specific tumors including brain, breast, colorectal, lung, thyroid, and others; please see separate listings for descriptions. The NeoTYPE Other Solid Tumor Profile is recommended for patients with resistant disease to explore options in clinical trials. This test is not for determining a tumor's tissue of origin.

### Specimen Requirements

- **FFPE tissue:** Paraffin block preferred. Please use 10% buffered formalin fixative. Do not use zinc fixatives.

### Storage & Transportation

Use cold pack for transport, making sure cold pack is not in direct contact with specimen.

### CPT Code(s)\*

81445x1, 81479x1, 88377x2, 88360x1, 88342x1; add 81479x1 if reflexed to NTRK NGS Fusion Panel (default) or 88374x3 automated (88377x3 manual) if reflexed to NTRK 1-3 FISH Panel

### Medicare MoIDX CPT Code(s)\*

81479

**New York Approved**

Yes

**Level of Service**

Global

**Turnaround Time**

14 days

NeoGenomics Laboratories is a specialized oncology reference laboratory providing the latest technologies, testing partnership opportunities, and interactive education to the oncology and pathology communities. We offer the complete spectrum of diagnostic services in molecular testing, FISH, cytogenetics, flow cytometry, and immunohistochemistry through our nation-wide network of CAP-accredited, CLIA-certified laboratories.

Committed to research as the means to improve patient care, we provide Pharma Services for pharmaceutical companies, in vitro diagnostic manufacturers, and academic scientist-clinicians. We promote joint publications with our client physicians. NeoGenomics welcomes your inquiries for collaborations. Please contact us for more information.

\*The CPT codes provided with our test descriptions are based on AMA guidelines and are for informational purposes only. Correct CPT coding is the sole responsibility of the billing party.

Please direct any questions regarding coding to the payor being billed.



9490 NeoGenomics Way  
Fort Myers, FL 33912  
Phone: 239.768.0600/ Fax: 239.690.4237  
neogenomics.com

© 2024 NeoGenomics Laboratories, Inc. All Rights Reserved.  
All other trademarks are the property of their respective owners  
Rev. 052024