

Nicolas M. Guibert^{1,2}, Yuebi Hu², Nora Feeney¹, Michael Epstein³, Samuel Woodhouse³, Karen Howarth³, Vincent Plagnol³, John Beeler³, Greg Jones³, Cloud P. Paweletz¹, Geoffrey R. Oxnard²

¹ Translational Research Laboratory; Belfer Center for Applied Cancer Science, Dana-Farber Cancer Institute; ² Lowe Center for Thoracic Oncology, Dana Farber Cancer Institute; ³ Inivata (Research Triangle Park, NC)


Background

- Several studies have evaluated hybrid-capture NGS for cfDNA genotyping, but this approach has some limitations¹
- Amplicon-based NGS is an attractive alternative with the potential to be faster and less expensive
- We performed a blinded evaluation of this approach for the characterization and monitoring of the molecular profile of advanced NSCLC during genotype-directed therapy

Methods

- 170 specimens from 47 pts were studied:
 - 145 from 31 advanced EGFR mutant NSCLC with T790M acquired resistance to EGFR-TKI
 - 25 specimens from 16 patients with other rare genotypes
- ddPCR for EGFR mutations was performed previously as part of ongoing research (DF/HCC protocol #14-147)²
- For each case, up to 8 specimens were analyzed, blinded to tumor genotype:
 - Pre-treatment
 - Initial follow-up draw
 - Subsequent follow-up draw(s)
 - Progression on TKI
- Plasma NGS was performed using the Inivata Liquid Biopsy Platform, Invision^{TM3,4}, which uses a combination of tagged amplicon sequencing and statistically-based analysis algorithms to identify and quantify low frequency tumor-derived single nucleotide variants (SNVs), short insertions/deletion (indels), copy number variants (CNVs) as well as structural variants of EML4/ALK fusions and ROS1 fusions. 36 cancer-related genes are examined using gene specific primers designed to hotspots and entire coding regions of interest.
- Plasma NGS results were reported blinded to any existing information on tumor and plasma genotyping
- Diagnostic accuracy was compared to tumor genotyping (including tumor NGS when available)

Inivata NGS panel



Genes established to be relevant to EGFR acquired resistance

EGFR	MET	BRAF	ERBB2	PIK3CA	PTEN
------	-----	------	-------	--------	------

Other key genes relevant to lung cancer biology

ALK	ROS1	RET	KRAS	MYC	FGFR1
FGFR2	AKT1	FGFR3	NRAS	TP53	STK11

Other cancer-related genes

GNAS	HRAS	CTNNB1	U2AF1	FOXL2	CCND1
IDH1	GNA11	GNAQ	KIT	MED12	ESR1
IDH2	PP2R1A	PDGFRA	NFE2L2	CDKN2A	NTRK3


Results

A. Amplicon-based Plasma NGS Is Highly Sensitive For The Detection Of A Full Range Of Genotypes

Tissue Genotype	Driver (AF%)		T790M(AF%)	
	ddPCR	NGS	ddPCR	NGS
	n=26		n=26	
L858R + T790M	1.1	0.65	0	0.28
L858R + T790M	0.7	1	0.45*	0
Del19 + T790M	0.65	0.47	0	0
L858R + T790M	0.5	0.23	0.4	0.57
L858R + T790M	0	0.057	0	0
Del19 + T790M	0	0.2	0.35*	0
Del19 + T790M	0	0.05	0	0
Sensitivity	88.5%	100%	84.5%	77%

Sensitivity of plasma NGS and plasma ddPCR for the detection of low AF EGFR mutations, compared to tissue. All other variants were detected by both ddPCR and NGS.

*low level positive, below threshold for clinical reporting



B. False Positives Are Rare With Amplicon-based Plasma NGS

Tables list all mutations detected in genes besides the driver, limited to genes regions sequenced with both assays

Driver	Baseline tissue NGS	Baseline plasma NGS
EGFR	TP53 A161T (55%)	none
EGFR	TP53 R110L (56%)	TP53 R110L (7.5%)
EGFR	TP53 166_167GA>A (64%)	TP53 frameshift (2.3%)
EGFR	none	none
EGFR	TP53 F134L (37%)	TP53 F134L (6.6%)
EGFR	none	PIK3CA E545K (0.6%)*
EGFR	TP53 L43* (32%)	TP53 L43 (8%)
EGFR	TP53 C242F (76%)	TP53 C242F (7.8%)
BRAF	STK11 p.D343N (30%)	STK11 p.D343N (50%)
MET	TP53 R248W	TP53 R248W (0.6%)
ROS1	none	none
ROS1	none	none
ALK	none	none
ALK	TP53 R337C	TP53 R337C (1.2%)




Tissue Genotype (SNVs)	NGS
BRAF V600E	KRAS G12D (0.12)**
BRAF V600E	none
BRAF G469A	BRAF G469A (12.8)
MET exon 14	MET exon 14 (1.4)
MET exon 14	MET exon 14 (12.1)
HER2 exon 20 ins	HER2 exon 20 ins (0.3)
Sensitivity	66.7%

Sensitivity for the detection of rare genotypes
** secondarily tested negative by ddPCR for KRAS and BRAF

A: Example of a CD74/ROS1 fusion gene detected in cfDNA

B, C: Examples of a EML4/ALK fusion gene detected in cfDNA

Plasma NGS can detect structural variants of EML4/ALK and ROS1 fusion genes with a 89% sensitivity (6/7 ALK, 2/2 ROS1).



Composite specificity calculated at 99.6% (i.e. the maximum number of false positive calls was 3 of 684 genes sequenced)

*The one suspected false positive, PIK3CA E545K, was then confirmed with ddPCR (AF 0.6%) and could be due to tumor heterogeneity

**Detection of a IDH1 mutation at low and constant AF, not following the same kinetics of other tumor-derived mutations, suggests a hematopoietic origin (e.g. CHIP).

C. Detection Of Acquired Resistance Mechanisms, Including Copy Number Variations

EGFR Patients with loss of T790M (n=7/15 (47%))


driver	Baseline NGS	Resistance NGS
Del19 Fig. A,D	BRAF V600E (0.60%)	BRAF V600E (12.4%) HER2 amp
Del19	BRAF V600E (6.5%) PTEN del (0,2%)	BRAF V600E (AF: 0.4%)
Del19 Fig. E	PIK3CA E545K (2.75%)	PIK3CA E545K (1.17%)
Del19 Fig. B	none	MET amp
Del19 Fig. C	PIK3CA E545K (0.6%)	HER2 amp
Del19	none	FGFR1 amp
L858R	none	KRAS G12S (0.21%)

EGFR Patients with maintained T790M (n=8/10, 80%)

Driver	Baseline NGS (beside driver and T790M)	Resistance NGS (beside driver and T790M)
Del19 Fig. F	none	EGFR C797S 26.6%
Del19 Fig. G	KRAS Q61K (0.12%)	EGFR C797S 7.9% KRAS Q61K (2.8%)
Del19	none	EGFR C797S 0.4% PTEN Y27C (21.9%) KRAS G13D (0.25%)
L858R	none	EGFR C797S_T_A 1.27% EGFR C797S_G_C 1.11%
L858R	none	EGFR C797S_T_A 2.7% EGFR C797S_G_C 1% EGFR Q791P_A_C 2.6%
del19	none	BRAF V600E 0.4% EGFR C797S_T_A 13.5%
del19	none	EGFR C797S_T_A 2.9%
del19	none	EGFR C797S_T_A 6.2% EGFR C797S_G_C 0.9%


Patients with rare genotypes (n=1/7 (14%))

Driver	Baseline NGS	Resistance NGS
ALK	EML4-ALK	ALK p.C1156Y 0.55



Patients with acquired amplifications at the time of resistance to osimertinib. A: Patient with acquired ERBB2 and MYC amplification (driver AF 63%)
B: Patient with acquired ERBB2 amplification and TP53 deletion at resistance (driver AF 38%). C: Patient with EGFR amplification and acquired MET amplification (driver AF 22%)

D. Early Detection Of Resistance in EGFR patients



Patients with mixed response, defined by rapid and complete response on the T790M clone and incomplete response on the driver. Plasma NGS can detect at baseline competing mutations: BRAF V600E mutation in patient D; PIK3CA E545K mutation in patient E.

Patients with acquired EGFR C797S tertiary mutations detected several months before progression. Plasma NGS additionally early detect a competing KRAS mutations in patient G. This KRAS mutation can be detected at low AF at baseline and re-emerge at resistance.

Conclusion

- In this blinded retrospective validation, amplicon-based plasma NGS could detect a full range of genotypes, including gene rearrangements with an exquisite sensitivity and a high specificity.
- This is in particular the first report of the detection of ALK and ROS1 fusions using amplicon sequencing, with a high sensitivity
- Quantitative concordance with ddPCR is excellent
- Serial plasma NGS can detect emergence of common resistance mutations
- Acquired KRAS mutation appears to be a newly identified recurring resistance mechanism to osimertinib
- In a subset of patients, resistance mutations can be detected at baseline, creating an opportunity for the study of targeted therapy combinations

Acknowledgments

This research was supported in part by the Expects Miracle Foundation, the Damon Runyon Cancer Research Foundation, the National Cancer Institute at the National Institutes of Health (R01 – CA135257), the U.S. Department of Defense (LCRP Career Development Award), the Anna Fuller Fund, and the Stading-Younger Cancer Research Foundation Thoracic Oncology Fund.

References

- Schwaederle MC, Patel SP, Husain H, Ikeda M, Lanman RB, Banks KC, et al. Utility of Genomic Assessment of Blood-Derived Circulating Tumor DNA (ctDNA) in Patients with Advanced Lung Adenocarcinoma. Clin Cancer Res. 2017 Sep 1;23(17):5101–11.
- Sacher et al. Prospective Validation of Rapid Plasma Genotyping for Detection of EGFR and KRAS Mutations in Advanced Lung Cancer. JAMA Oncol. 2016;2:1014–22.
- Forshew T, Murtaza M, Parkinson C, Gale D, Tsui DW, Kaper F, et al. Noninvasive identification and monitoring of cancer mutations by targeted deep sequencing of plasma DNA. Sci Transl Med. 2012 May 30;4(136):136ra68
- Gale D, Plagnol V, Lawson A, Pugh M, Smalley S, Howarth K, et al. Cancer Research, AACR 2016, Abstract 3639: Analytical performance and validation of an enhanced TAM-Seq circulating tumor DNA sequencing assay; DOI: 10.1158/1538-7445.AM2016-3639 Published July 2016