



# Test Catalog

Diagnostic. Prognostic. Predictive. Predisposition.



# Amyloid A & Amyloid P

## Alternative Name

Amyloid A, Amyloid P, Amyloid A&P Panel

## Methodology

Immunohistochemistry (IHC)

## Test Description

Amyloid A and P react with amyloid deposits in many tissues. When accompanied by Congo Red, Amyloid A and P can be used to distinguish primary and secondary amyloidosis. Because these stains are interpreted in context of each other and Congo Red, testing options available to global and tech-only clients differ.

## Specimen Requirements

- **Global orders:** Please order Amyloid A&P Panel and submit a Congo Red-stained slide for our reference (required). If Congo Red staining was not done or the slide is not available, please order a Consult instead so we can coordinate Congo Red staining followed by amyloid A and P staining as indicated. Individual orders for Amyloid A or Amyloid P are not available with global service.
- **Tech-only (stain only) orders:** Please order Amyloid A and/or Amyloid P as individual stains. The Panel combination is not available with tech-only service, nor is it necessary to submit a Congo Red-stained slide.
- A formalin-fixed, paraffin-embedded (FFPE) tissue block is preferred specimen type  
or
- One (1) unbaked, unstained slide for H&E staining (required) and two to four (2-4) positively charged unstained slides (all cut at 4-5 microns)
- Block and slide identifiers should be clearly written and match exactly with the specimen ID and specimen labeling as noted on the requisition.

## Storage & Transportation

Use cold pack for transport, making sure cold pack is not in direct contact with specimen.

## CPT Code(s)\*

Panel for global clients: (88342x1, 88341x1) or 88341x2. Single stain for tech-only clients: 88342x1 or 88341x1 for each stain

## New York Approved

Yes

## Level of Service

Global, Stain Only

## Turnaround Time

Global: 48 hours, Tech-Only (stain only): 24 hours

\*The CPT codes provided with our test descriptions are based on AMA guidelines and are for informational purposes only. Correct CPT coding is the sole responsibility of the billing party.

Please direct any questions regarding coding to the payor being billed.

NeoGenomics Laboratories is a specialized oncology reference laboratory providing the latest technologies, testing partnership opportunities, and interactive education to the oncology and pathology communities. We offer the complete spectrum of diagnostic services in molecular testing, FISH, cytogenetics, flow cytometry, and immunohistochemistry through our nation-wide network of CAP-accredited, CLIA-certified laboratories.

Committed to research as the means to improve patient care, we provide Pharma Services for pharmaceutical companies, in vitro diagnostic manufacturers, and academic scientist-clinicians. We promote joint publications with our client physicians. NeoGenomics welcomes your inquiries for collaborations. Please contact us for more information.

\*The CPT codes provided with our test descriptions are based on AMA guidelines and are for informational purposes only. Correct CPT coding is the sole responsibility of the billing party.

Please direct any questions regarding coding to the payor being billed.



9490 NeoGenomics Way  
Fort Myers, FL 33912  
Phone: 239.768.0600/ Fax: 239.690.4237  
neogenomics.com

© 2024 NeoGenomics Laboratories, Inc. All Rights Reserved.  
All other trademarks are the property of their respective owners  
Rev. 052024