



Test Catalog

Diagnostic. Prognostic. Predictive. Predisposition.





Albumin RNA ISH

Methodology

In Situ Hybridization (ISH)

Test Description

Albumin RNA ISH is a sensitive and specific tool for distinguishing primary hepatocellular carcinoma and intrahepatic cholangiocarcinoma from metastatic adenocarcinoma to the liver or carcinoma of unknown origin in formalin fixed paraffin-embedded tissues. Positive results in this assay provide evidence of liver origin.

Specimen Requirements

- A formalin-fixed, paraffin-embedded (FFPE) tissue block is preferred specimen type
or
- One (1) unbaked, unstained slide cut at 4-5 microns for H&E staining (required) and three (3) positively charged unstained slides, all cut at 4-5 microns
- Block and slide identifiers should be clearly written and match exactly with the specimen ID and specimen labeling as noted on the requisition

Storage & Transportation

Use cold pack for transport. Make sure cold pack is not in direct contact with specimen.

CPT Code(s)*

88365x1

New York Approved

Yes

Level of Service

Stain Only

Turnaround Time

48 hours

*The CPT codes provided with our test descriptions are based on AMA guidelines and are for informational purposes only. Correct CPT coding is the sole responsibility of the billing party.

Please direct any questions regarding coding to the payor being billed.

NeoGenomics Laboratories is a specialized oncology reference laboratory providing the latest technologies, testing partnership opportunities, and interactive education to the oncology and pathology communities. We offer the complete spectrum of diagnostic services in molecular testing, FISH, cytogenetics, flow cytometry, and immunohistochemistry through our nation-wide network of CAP-accredited, CLIA-certified laboratories.

Committed to research as the means to improve patient care, we provide Pharma Services for pharmaceutical companies, in vitro diagnostic manufacturers, and academic scientist-clinicians. We promote joint publications with our client physicians. NeoGenomics welcomes your inquiries for collaborations. Please contact us for more information.

*The CPT codes provided with our test descriptions are based on AMA guidelines and are for informational purposes only. Correct CPT coding is the sole responsibility of the billing party.

Please direct any questions regarding coding to the payor being billed.



9490 NeoGenomics Way
Fort Myers, FL 33912
Phone: 239.768.0600/ Fax: 239.690.4237
neogenomics.com

© 2024 NeoGenomics Laboratories, Inc. All Rights Reserved.
All other trademarks are the property of their respective owners
Rev. 052024